

# CAAT

CARBONACEOUS AEROSOL • IMPORT your data from Magee • Analyze and display the data ANALYSIS TOOL

Advanced Black Carbon Measurement

## **KEY FEATURES**

- Scientific instrument
- VALIDATE your data.
- ANALYZE your validated data by using several procedures: Diurnal profile, Time Average, Absorption Coefficient, Ångstrom exponent, Attenuation, Source Apportionment.
- DISPLAY your data by using various chart options.

## **APPLICATIONS**

from Aerosol Magee Scientific instruments

# AEROSOL MAGEE SCIENTIFIC

#### PRODUCT SPECIFICATIONS

Analyze your Magee Scientific instrument data with CAAT Software.

#### BENEFITS

Be assured that your data is correct and consistent Use CAAT as your primary tool for your analysis of data obtained from Magee Scientific instruments. Immediately check the proper "data status". Start using CAAT Software tool No additional requirements. Straightforward data evaluation and display. User friendly controls.


#### **EASY INSTALLATION**

- Login or Register at Magee Scientific User's group: Search for CAAT product.
- Go to SOFTWARE tab and download installation files.
- Search for the downloaded files in your computer.
  Unzip to a selected destination folder.
- Click on the "CAAT-Setup.exe" and follow the guided procedure.
- You can start using the CAAT application immediately after installation is completed.

#### CAAT







#### AEROSOLMAGEESCI.COM





Address: U5 / 17 Enterprise Crescent, McDougalls Hill NSW 2330 Australia